

## SRT-Tag Mouse Monoclonal Antibody(4H9)

Catalog No: #T518

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SRT-Tag Mouse Monoclonal Antibody(4H9)
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	N/A
Specificity	The SRT Tag Mouse Monoclonal Antibody detects SRT Tag recominant protein.
Target Name	SRT-Tag
Other Names	SRT Tag antibody
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 51% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:5000~1:10000

## Images



Western blot analysis of SRT-Recombinant protein, using #T518 diluted at 1) 1:5,000 2) 1:10,000

## Background

The SRT-tag is a peptide epitope comprised of 10 amino acids(TFIGAIATDT) derived from the crystalline surface layer protein of Rickettsia typhi1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.