

T7-Tag Mouse Monoclonal Antibody

Catalog No: #T604

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	T7-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB IF/ICC
Species Reactivity	Tag
Specificity	T7-tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal T7-tagged proteins.
Immunogen Type	Peptide
Immunogen Description	A synthetic peptide MASMTGGQQMG
Target Name	T7-Tag
SDS-PAGE MW	N/A
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C. Stable for 12 months from date of receipt

Application Details

WB: 1:3000~1:10000

IF/ICC: 1:200

Images



Western blot analysis of GST-T7 fusion protein, using T7-tag mouse monoclonal antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.