

## RFP-Tag Mouse Monoclonal Antibody

Catalog No: #T608

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RFP-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB IHC IF/ICC
Species Reactivity	Tag
Specificity	RFP-Tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal RFP recombinant proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Full length RFP protein.
Target Name	RFP-Tag
SDS-PAGE MW	N/A
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C. Stable for 12 months from date of receipt


## Application Details

WB: 1:3000~1:10000

IHC: 1:200

IF/ICC: 1:200

## Images



1ug RFP fusion protein+ Primary antibody dilution at  
 1γ • 1:3,000 2γ • 1:5,000 3γ • 1:10,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.