

Myc-Tag Monoclonal Antibody

Catalog No: #T612



Package Size: #T612-1 50ul #T612-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Myc-Tag Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	N/A
Specificity	The antibody can detect C-Myc fusion proteins.
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:5000-10000

Images



Western blot analysis of His-MBP-C-Myc recombinant protein, diluted at
1) 1:5000
2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.