
HbA1C - low level in normal Hb Conc. Datasheet

Product details

- **Product Code:** PS186-0
- **Type:** Native Antigens
- **Applications:** Control Manufacture Life Science Clinical Chemistry Biosensors ELISA Assay Lateral Flow
- **Shipping recommendation:** 4°C

Product summary

- **Presentation matrix:** stabiliser
- **Source:** human blood
- **Abbreviation:** HbA1C low

Application notes

- **Recommended Use:** Hba1c controls solns

Additional Information

- **Shipping recommendation:** 4°C
- **Storage:** 4°C if thawed
- **Keywords:** diabetic, HbA1C, diabetes, hemoglobin, haemoglobin, A1C