



Certificate of Analysis

Product: **Goat anti-Mouse IgG (H+L)-PolyHRP40, polyclonal affinity purified Goat IgG antibody, specific to heavy and light chains of mouse IgG, conjugated with PolyHRP40**

Code #: **#GAMP40C**

Concentration: **1 mg/ml total conjugate**

Buffer base: **0,1M K-Phosphate buffer pH 6,0 + 50% ultra-pure glycerol**

Stabilizing additives: HRP stabilizers (proprietary) - 100 & 250 µg/ml  
Ultra pure IgG-free BSA - 15 mg/ml

Preservative/anti-microbial: 5-Bromo-5-nitro-1,3-dioxane - 600 ppm

Lot: 222638

Storage: -20°C (-15°C/-25°C). Do not store below -25°C!

Expiration date: 10/2023

Quality Control Assurance Release Procedure: ELISA use test with rabbit IgG as an antigen

Result: **POSITIVE (COMPLIES)** ~~NEGATIVE (DOES NOT COMPLY)~~

Date: 23.09.2020