

Certificate of Analysis

Product: **Goat anti-Mouse IgG (H+L)-PolyHRP40,**

Code number: **#GAMP40D50**

Concentration: **50 µg/ml pre-diluted total conjugate in liquid formulation stabilized for storage at +4°C (format 1/20)**

Buffer base: **Antibody/Antigen (PolyHRP) Stabilizer #AA1**

Stabilizing additives: **proprietary**

Preservative/anti-microbial: **5-Bromo-5-Nitro-1,3-Dioxane - 1200 ppm**

Lot: **766129**

Storage: **+4°C (+2°C/+8°C/+25°C). Do not freeze!**

Expiration date: **05/2014**

Quality Control Assurance Release Procedure: **ELISA use test with maus IgG as an antigen at different conditions**

Result:

**POSITIVE
(COMPLIES)**

~~**NEGATIVE
(DOES NOT COMPLY)**~~

Date:

25.04.2013

This certificate is a computer printout and has therefore not been signed by hand