

Manufacturer of Specialty Reagents for ImmunoAssays, DNA/RNA-hybridization and ImmunoHistoChemistry

Arnold-Sommerfeld-Ring 2 52499 Baesweiler • Germany Tel.: 02401-606870

Fax: 02401-608534 support@sdt-reagents.de

## **Certificate of Analysis**

| Product:  | Goat anti-Rabbit IgG (H+L)-PolyHRP40  |
|---|---|
| Code#:  | #GARP40D50  |
| Concentration:                                  | 50 $\mu$ g/ml pre-diluted <b>total conjugate</b> in liquid formulation stabilized for storage at +4°C (format 1/20) |
| Buffer base:                                    | SA-PolyHRP Stabilizer #SA1  |
| Lot:  | 416463  |
| Storage:  | $+2^{\circ}C/+8^{\circ}C$ . Do not freeze below $0^{\circ}C!$   |
| Expiration date:                                | 05/2014   |
| Quality Control Assurance<br>Release Procedure: | ELISA use test with rabbit IgG as analyte at different conditions   |
| Result:   | POSITIVE<br>(COMPLIES) NEGATIVE<br>(DOES NOT COMPLY)  |
| Date:   | 25.04.2013  |

This certificate is a computer printout and has therefore not been signed by hand