



Certificate of Analysis

Product: **Streptavidin-PolyHRP40-DY547**

Code #: **#SP40-DY547C**

Concentration: **1 mg/ml total conjugate**

Buffer base: **PBS pH 7,3 + 50% ultra-pure glycerol**

Stabilizing additives: HRP stabilizer (proprietary) - 100 µg/ml
 Ultra pure biotin-free BSA - 15 mg/ml

Preservative/anti-microbial: 5-Bromo-5-nitro-1,3-dioxane - 600 ppm

Lot: **386733**

Storage: **-20°C (-15°C/-25°C). Do not store below -25°C!**

Expiration date: **04/2023**

Quality Control Assurance ELISA use test with hIgG-biotin as analyte
Release Procedure: at different conditions

Result: **POSITIVE**
 (COMPLIES)

NEGATIVE
 (DOES NOT COMPLY)

Date: **19.03.2020**