



Certificate of Analysis

Product: **Streptavidin-PolyHRP80**

Code #: **#SP80C**

Concentration: **1 mg/ml total conjugate**

Buffer base: **PBS pH 7,3 + 50% ultra-pure glycerol**

Stabilizing additives: **HRP stabilizer (proprietary) - 100 µg/ml**
Ultra pure biotin-free BSA - 15 mg/ml

Preservative/anti-microbial: **5-Bromo-5-nitro-1,3-dioxane - 600 ppm**

Lot: **417983**

Storage: **-20°C (-15°C/-25°C). Do not store below -25°C!**

Expiration date: **12/2024**

Quality Control Assurance
Release Procedure: **ELISA use test with hIgG-biotin as analyte
at different conditions**

Result: **POSITIVE
(COMPLIES)** ~~**NEGATIVE
(DOES NOT COMPLY)**~~

Date: **11.11.2021**