

Manufacturer of Specialty Reagents for ImmunoAssays,  
DNA/RNA-hybridization and ImmunoHistoChemistry

Fax: 02401-608534  
support@sdt-reagents.de

**Certificate of Analysis**

Product: **Streptavidin-PolyHRP80**

Code #: **#SP80D50**

Concentration: **50 µg/ml pre-diluted SA-PolyHRP80, liquid formulation stabilized for storage at +4°C (format 1/20)**

Buffer base: **Streptavidin-PolyHRP Stabilizer #SA1-HTS**

Lot: **765727**

Storage: **+2°C/+8°C. Do not freeze below 0°C!**

Expiration date: **11/2014**

Quality Control Assurance  
Release Procedure: **ELISA use test with hIgG-biotin as analyte at different conditions**


Result:

**POSITIVE  
(COMPLIES)**

~~**NEGATIVE  
(DOES NOT COMPLY)**~~

Date:

**06.11.2013**



Stereospecific  
Detection Technologies

Arnold-Sommerfeld-Ring 2  
52499 Baesweiler • Germany  
Tel.: 02401-606870

---

Manufacturer of Specialty Reagents for ImmunoAssays,  
DNA/RNA-hybridization and ImmunoHistoChemistry

Fax: 02401-608534  
support@sdt-reagents.de

This document has been electronically produced and is valid without a signature