


### pAcYM1 Baculovirus Transfer Vector

#### Product Summary

**Cat. No:** STV01  
**Quantity:** 10 µg/20 µl  
**Storage:** Store at -20°C  
**Stability:** Stable for 6 months in appropriate conditions

#### Description:

pAcYM1 is a baculovirus transfer vector designed for high level expression of foreign genes under the powerful AcMNPV polyhedrin promoter. The transfer vector is the modified transfer vector and is used to insert the foreign gene in the cloning steps. It has a bacterial origin of replication and an ampicillin resistance gene in pUC8. The polyhedrin gene sequences are deleted and there are both 5' noncoding and 3' noncoding DNA sequences of polyhedrin gene in pAcYM1. BamHI is used for insertion of the foreign gene and restriction sites, including XhoI, EcoRV and HindIII, are used to confirm the orientation. The AcMNPV sequences flanking the gene, 4kb and 2kb, in the transfer vectors allow homologous recombination with the viral DNA to insert the expression cassette into the polyhedrin locus. pAcYM1 is compatible with superBAC derivatives and any baculovirus system that utilizes homologous recombination in sf9 insect cells.



### pAcYM1 polyhedrin promoter sequence

GATATCATGGAGATAATTTAAAATGATAACCATCTCGCAAATAAA  
EcoRV

TAACTATTTTACTGTTTTTCGTAACAGTTTTGTAATAAAAAAACCC  
polyhedrin mRNA 5' end

TATAAATCGGATCCGTTATT  
BamHI

### Related Products

Product	Cat. No.
pAcYM1_L21	STV02
pAc	STV03
pAc_L21	STV04
pAcYM1-polyhedrin baculovirus control vector	STV11
Baculovirus genomic DNA with polyhedrin gene	STV21
superBAC_wt baculovirus DNA	ST001
superBAC_PDI baculovirus DNA	ST002
superBAC_cdc37 baculovirus DNA	ST003

**Please contact us for a complete list of the baculovirus related research products, services and further information.**

#### SHEATECH, INC.

6339 Nancy Ridge Drive, Suite 200  
 San Diego, CA 92121  
 Tel: (619) 788-8188

Email: [info@sheatechinc.com](mailto:info@sheatechinc.com)  
<http://www.sheatechinc.com>

**For Research Use Only. Not for Diagnostic or Therapeutic Use.**

This transfer vector is used for cloning the foreign gene only. It is not used for diagnostic or therapeutic use. Any use of this product other than the permitted use without the express written authorization of SHEATECH is prohibited.