


Product Information

Product Name	PVLA
Structure	
Product Molecular Mass	Mw:170,000
Lot No.	44167
Reagent	powder
Concentration	—
Quantity	100mg
Stability / Storage	Store at room temperature.
Purity	85%
Protocol	Prepared at time of use 1. Dissolve PVLA in water to 100 µg / ml. 2. Pass through a sterilization filter. 3. Coat the culture vessel and incubate at 37°C for 2hr. 4. Washed with PBS three times, and then cells were seeded.

For research use only