

TRAIL Antibody : B-S23 Azide free

ANTI HUMAN TRAIL/CD253 AZIDE FREE

CATALOGUE N° : 853.080.000

CLONE : B-S23

ISOTYPE : Mouse IgG1

IMMUNISATION : Recombinant human TRAIL

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the TRAIL, APO2L antigen, a 32 kDa protein

WORKSHOP : VIII

QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Suitable for flow cytometry
Enhances TRAIL induced apoptosis
Immunoprecipitation : not tested
Western Blot : not tested
Immunohistochemistry : not tested
Recognises the soluble antigen

REFERENCES : Vermot-Desroches C. et al., Cell Immunol.,
2005;236(1-2):86-91