

CD117 Antibody : B-K15 Azide free

HUMAN CD117 AZIDE FREE

CATALOGUE N° : 854.510.000

CLONE : B-K15

ISOTYPE : Mouse IgG1

IMMUNISATION : Recombinant c-kit

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the Stem Cell Factor Receptor, SCF R, c-kit, a 145 kDa protein


QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months. For longer storage, freeze aliquots.

APPLICATION : Suitable for flow cytometry
Inhibits the Kit-Ligand induced proliferation of TF1 cell line
Immunoprecipitation : not tested
Western Blot : not tested
Immunohistochemistry : not tested
Recognises the soluble antigen



A typical staining pattern with the B-K15 monoclonal antibody of TF1 cell line