

CD184 Antibody : B-R24 Azide free

HUMAN CD184 AZIDE FREE

CATALOGUE N° : 854.600.000

CLONE : B-R24

ISOTYPE : Mouse IgG1

IMMUNISATION : CXCR4 transfected CHO cell line

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the CXCR4 antigen, a 39.7 kDa protein

WORKSHOP : VII


QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months. For longer storage, freeze aliquots

APPLICATION : Suitable for flow cytometry
Inhibits HIV infection
Immunoprecipitation : not tested
Western Blot : not tested
Immunohistochemistry : not tested



A typical staining pattern with the B-R24 monoclonal antibody of lymphocytes