

CD202b Antibody : B-B48 Azide free

HUMAN CD202b AZIDE FREE

CATALOGUE N° : 854.710.000

CLONE : B-B48

ISOTYPE : Mouse IgG1

IMMUNISATION : Eahy 926 cell line

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises Tie2, Tek, an unknown KDa protein

WORKSHOP : VII


QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months. DO NOT FREEZE

APPLICATION : Suitable for flow cytometry - Inhibits Tek-L / Tek binding
Immunoprecipitation : not tested - Western Blot : not tested
Immunohistochemistry : not tested



A typical staining pattern with the B-B48 monoclonal antibody of Eahy 926 cell line