



## TRAIL R3 Antibody : B-D44 Azide free

ANTI HUMAN TRAIL R3/CD263 AZIDE FREE

**CATALOGUE N° :** 854.890.000

**CLONE :** B-D44

**ISOTYPE :** IgG1

**IMMUNISATION :** Recombinant human TRAIL R3/Fc chimera

**HYBRIDOMA :** Myeloma X63/AG.8653 x Balb/c spleen cells

**SPECIFICITY :** Recognises the 75.90kDa protein

**QUANTITY :** 200 µg / 200 µl

**PURIFICATION :** Ion exchange chromatography

**FORMAT :** Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

**STORAGE :** Stable at +2-8°C for 12 months  
For longer storage, freeze aliquots

**APPLICATION :** Inhibits TRAIL binding to TRAIL R3

**REFERENCES :** Vermot-Desroches C. et al., Cell Immunol.,  
2005;236(1-2):86-91

**Version 3 - 01.08**

**For research use only**