

CDw218a Antibody : B-E43 Azide free

HUMAN CDw218a/IL-18Ra AZIDE FREE

CATALOGUE N° : 854.900.000

CLONE : B-E43

ISOTYPE : Mouse IgG1

IMMUNISATION : Recombinant human IL-18Ra/Fc

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the IL-18Ra, a 63 kDa protein

WORKSHOP : VIII

QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Blocks human IFNγ secretion induced by IL-18 stimulated cells

REFERENCES : Vermot-Desroches C. et al., Cell Immunol., 2005; 236 (1-2):101-104

Version 2 - 01.06

For research use only