

CDw218b Antibody : B-B46 Azide free

HUMAN CDw218b/IL-18Rb AZIDE FREE

CATALOGUE N° : 854.910.000

CLONE : B-B46

ISOTYPE : Mouse IgG1

IMMUNISATION : Recombinant human IL-18Rb/Fc

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the IL-18Rb, a 65.1 kDa protein

QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Blocks human IFNγ secretion induced by IL-18 stimulated cells

REFERENCES : Vermot-Desroches C. et al., Cell Immunol., 2005;
236(1-2):101-104

Version 2 - 01.06

For research use only