

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Monkey IL-1 β ELISA kit

Catalogue no.: CT139A (5-plate format)

Validated for:

Rhesus macaque and Barbary macaque



Recommended assay range: 8-256 pg/ml.

Sensitivity: 5 pg/ml.

Specificity: recombinant and natural rhesus macaque IL-1 β .

Calibration: recombinant rhesus macaque IL-1 β .

Type of sample: serum, plasma and culture supernatant.