


Typical data Monkey IL-6 ELISA kit

Catalogue no.: CT145A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Barbary macaque, Japanese macaque, Baboon, African green monkey



Recommended assay range: 1-32 pg/ml.

Sensitivity: 1 pg/ml.

Calibration: natural rhesus macaque IL-6.

Type of sample: serum, plasma and culture supernatant.

Minor interference (<10%) by the sIL-6 receptor at a 1000-fold molar excess.

References list:

Beignon, A.-S., Mollier, K., Liard, C., Coutant, F., Munier, S., Rivière, J., Souque, P. and Charneau, P. 2009. *J. Virol.* 83: 10963-10974

Coe, C.L., Lubach, G.R. and Kinnard, J. 2012. *Age (Dordr)* 34: 1169-1177

McFarlane, D., Wolf, R.F., McDaniel, K.A. and White, G.L. 2012. *J. Med. Primatol.* 41: 341-348

Rout, N., Greene, J., Yue, S., O'Connor, D., Johnson, P.R., Else, J.G., Exley, M.A. and Kaur, A. 2012. *PloS Pathog.* 8: e1002928

More references can be found in our References Database on our website: www.ucytech.com/references.