

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Monkey IL-13 ELISA kit

Catalogue no.: CT147A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Japanese macaque, Lion-tailed macaque, Barbary macaque, Baboon, African green monkey, Marmoset



Recommended assay range: 0.5-32 pg/ml.

Sensitivity: 0.5 pg/ml.

Calibration: natural rhesus macaque IL-13.

Type of sample: serum, plasma and culture supernatant.

References list:

Rout, N., Greene, J., Yue, S., O'Conner, D., Johnson P.R., Else, J.G., Exley, M.A. and Kaur, A. 2012. *PLoS Pathog.* 8: e1002928

De Swart, R.L., Kuiken, T., Timmermans, H.H., van Amerongen, G., van den Hoogen, B.G., Vos H.W., Neijens, H.J., Andeweg, A.C. and Osterhaus, A.D.M.E. 2002. *J.Virol.* 76: 11561-11569

De Waal, L., Power, U.F., Yüksel, S., van Amerongen, G., Nguyen, T.N., Niesters, H.G.M., de Swart, R.L. and Osterhaus, A.D.M.E. 2004. *Vaccine* 22: 915-922

Wang, C.Y., Finstad, C.L., Walfield, A.M., Sia, C., Sokoll, K.K., Chang, T-Y., De Fang, X., Hung, C.H., Hutter-Paier, B. and Windisch, M. 2007. *Vaccine* 25: 3041-3052

More references can be found in our References Database on our website: www.ucytech.com/references.