

U-CyTech BV


Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Monkey Perforin ELISA kit

Catalogue no.: CT154A (5-plate format)

### Validated for:

Rhesus macaque and Barbary macaque



Recommended assay range: 16-1024 pg/ml

Sensitivity: 12 pg/ml.

Calibration: natural rhesus macaque Perforin.

Type of sample: serum, plasma and culture supernatant.