


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Marmoset IFN- γ ELISA kit

Catalogue no.: CT340A (5-plate format)

Validated for:

Marmoset



Recommended assay range: 0.5-32 pg/ml.

Sensitivity: 0.5 pg/ml.

Calibration: recombinant and natural marmoset IFN- γ .

Type of sample: culture supernatant.

References list:

Jagessar, A.S., Heijmans, N., Bauer, J., Blezer, E.L.A., Laman, J.D., Migone, T.-S., Devalaraja, M.N. and 't Hart, B.A. 2012. *J. Neuroimmune Pharmacol.* 7: 557-570