


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Marmoset IL-13 ELISA kit

Catalogue no.: CT341A (5-plate format)

Validated for:

Marmoset



Recommended assay range: 0.5-16 pg/ml.

Sensitivity: 0.5 pg/ml.

Calibration: natural marmoset IL-13.

Type of sample: serum, plasma and culture supernatant.