


Typical data Monkey IL-17A ELISA kit

Catalogue no.: CT501A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey and Baboon



Recommended assay range: 0.5-64 pg/ml.

Sensitivity: 0.5 pg/ml.

Calibration: natural rhesus macaque IL-17A.

Type of sample: serum, plasma and culture supernatant.

References list:

Rout, N., Greene, J., Yue, S., O'Conner, D., Johnson P.R., Else, J.G., Exley, M.A. and Kaur, A. 2012. PLoS Pathog. 8: e1002928