

U-CyTech BV

Yalelaan 48  
3584 CM Utrecht  
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

## Typical data Monkey IL-23 ELISA kit

Catalogue no.: CT502A (5-plate format)

### Validated for:

Rhesus macaque, Cynomolgus monkey



Example of a representative Monkey IL-23 standard curve using a 4-parameter logistic fit.

Recommended assay range: 8-512 pg/ml.

Sensitivity: 8 pg/ml.

Calibration: natural rhesus macaque IL-23.

Type of sample: culture supernatant.