

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Monkey IL-17F ELISA kit

Catalogue no.: CT503A (5-plate format)

Validated for:

Rhesus macaque



Example of a representative Monkey IL-17F standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.

Sensitivity: 4 pg/ml.

Calibration: natural rhesus macaque IL-17F.

Type of sample: serum, plasma and culture supernatant.