


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Marmoset IL-6 CytoMax ELISA kit

Catalogue no.: CT974-2 (2-plate format)
CT974-5 (5-plate format)

Validated for:

Marmoset



Recommended assay range: 2.5 - 80 pg/ml.
Sensitivity: 2.5 pg/ml.
Specificity: natural marmoset IL-6.
Calibration: natural marmoset IL-6.
Type of sample: serum, plasma and culture supernatant.