



CERTIFICATE OF ANALYSIS

Product ALKALINE PHOSPHATASE AAL (ALEURIA AURANTIA LECTIN)
AFFINITY-PURIFIED

Catalog No. MB-4100

Amount 150 µg

Lot No. W0927

Concentration 0.5 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.
Immunohistochemical grade bovine serum albumin has been added
for stabilization.

Storage Refrigerate; DO NOT FREEZE

Usage Information:

For DNA/RNA transfer blot applications, A dilution of 1:5000 made in 1X Casein Solution (Cat. No. SP-5020) or TBS or PBS containing 0.1% Tween 20 is recommended with a suggested incubation time of 30 minutes at room temperature with BCIP/NBT as the substrate.

Please note that the above conditions and accompanying protocols are provided only as guidelines-conditions should be optimized for individual cases.

see reverse for additional applications