

Visikol® HISTO™ 10X Washing Buffer

Antigen Retrieval Media

INTRODUCTION

Visikol HISTO 10X Washing Buffer is an antigen retrieval buffer designed to aid in visualization of tissues in 3D. Using the Visikol HISTO 10X Washing Buffer gives researchers more time to focus on what matters, instead of preparing media.

TYPICAL APPLICATIONS

- Antigen retrieval agent in Visikol® HISTO™ tissue clearing protocols

ADVANTAGES

- Reduces reagent preparation time in tissue clearing protocols
- Ready to use

TYPICAL PHYSICAL PROPERTIES

These properties are typical but should not be considered specifications. Please contact us to obtain the certificate of analysis for your specific lot of material.

Appearance:	Clear, colorless liquid
	PBS
Components:	Tween® 20
	Heparin

DISCLAIMER

No freedom from infringement of any patent owned by Visikol or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other governmental enactments. The product shown in this literature may not be available for sale and/or available in all geographies where Visikol is represented. The claims made may not have been approved for use in all countries. Visikol assumes no obligation or liability for the information in this document. NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED.