

Mature BDNF Rapid™ ELISA Kit

Mature BDNF Rapid™ ELISA Kit (BEK-2211)

- Measure BDNF in serum, plasma, cell culture supernatants, cell lysates and brain extracts **in less than 3 hours!**
- For the preferential detection of **Mature BDNF** isoform independent of the proBDNF isoform.
- Quantify BDNF in **Human and Rodent** samples against a **WHO approved BDNF standard**.
- Kits contain a positive **Quality Control** sample & strip wells for ease of use and quality assurance.
- High sensitivity gives you the advantage of using small **sample volumes**.
- **Focus** on your experiment - Let us take the worry out of reliable quantification of BDNF concentrations.

“The Best BDNF kit” as
Published in...
Polacchini et al. (2015) Nature Sci.
Rep. 10;5:17989



**CE Marked
Version
Available!**

Cat# BEK-2211-1P/2P-CE

What forms of BDNF are detected:


Mature BDNF Rapid ELISA Kit Catalog No: BEK-2211

For the preferential detection of Mature BDNF independently of the proBDNF form.



Below figure: Our Mature BDNF Rapid™ ELISA Kit Cat# BEK-2211-1P/2P is specific for BDNF and shows little to no value when trying to measure proBDNF.

BEK-2211 Mature BDNF ELISA



Biosensis Mature BDNF Rapid™ ELISA

Target	Code (2 Plate)	Code (1 Plate)	Species reactivity	Sensitivity	Range
BDNF Rapid™	BEK-2211-2P	BEK-2211-1P	Hu, Ms, Rat	1-2 pg/mL	7.8 pg/mL - 500 pg/mL