



Start of Transcription

Polyhedrin promoter

4481 ATGGAGATAA TTAAAATGAT AACCATCTCG CAAATAAATA AGTATTTTAC
wild-type ATG mutated to ATT

4531 TGTTTTTCGTA ACAGTTTTGT AATAAAAAAA CCTATAAATA TTCCGGATTA

4581 TTCATACCGT CCCACCATCG GCGCGGATC CCGGTCCGAA GCGCGCGGAA
Bam H I Rsr II BssH II EcoR I

4631 TTCAAAGGCC TACGTCGACG AGCTCACTAG TCGCGGCCGC TTTCGAATCT
Stu I Sal I Sst I Spe I Not I Nsp V Xba I

4681 AGAGCCTGCA GTCTCGACAA GCTTGTCGAG AAGTACTAGA GGATCATAAT
Pst I Xho I Hind III

SV40 polyadenylation signal

4731 CAGCCATACC ACATTTGTAG AGGTTTACT TGCTTTAAAA AACCTCCCAC