



New Released Antibodies for Cell Signal transduction Research


Cat. No.	Product name	Reactivity	Applications
11553	Cdc25B (Phospho-Ser149) Antibody	Mouse Rat	WB
21437	CDK4 Antibody	Human	WB
21667	C-Fos Antibody	Human Mouse Rat	WB
11572	CK2 α (phospho-Thr360/Ser362) Antibody	Human, Mouse, Rat	WB
21572	CK2 α Antibody	Human, Mouse, Rat	WB
11240	c-Kit (Phospho-Tyr721) Antibody	Human Mouse Rat	WB
21679	Cox2 Antibody	Human	WB
21513	eIF4B (Ab-422) Antibody	Human	WB IHC
11233	eIF4E (Phospho-Ser209) Antibody	Human	WB IHC IF
21514	eIF4G (Ab-1232) Antibody	Human Mouse Rat	WB IHC IF
21619	Glut4 Antibody	Human Mouse Rat	WB
21613	Glycogen Synthase Antibody	Human Mouse Rat	WB
21660	HDAC3 Antibody	Human Mouse Rat	WB
11577	Histone H3 (Phospho-Thr11) Antibody	Human Mouse Rat	WB
11554	IL-10R subunit alpha (Phospho-Ser319/323) Antibody	Human	WB
11558	Interferon- α Receptor Type I Subunit I (phospho-Ser535/Ser539) Antibody	Human	WB
21681	IRF-3 (Ab-396) Antibody	Human	WB
21615	Keratin 19 Antibody	Human	WB
21441	MAP2K4 Antibody	Human Mouse Rat	WB
11559	MEK5(phospho-Ser311/ Thr315) Antibody	Human Mouse Rat	WB
11227	Met (Phospho-Tyr1234) Antibody	Human Mouse Rat	WB
11238	Met (Phospho-Tyr1349) Antibody	Human Mouse Rat	WB
21515	mTOR (Ab-2481) Antibody	Human Mouse Rat	IHC
11221	mTOR (Phospho-Ser2448) Antibody	Human	WB IHC IF
11578	Niban (Phospho-Ser602) Antibody	Human	WB
21683	P38 MAPK(Ab-180/182) Antibody	Human Mouse Rat	WB
11581	P38 MAPK(Phospho-Thr180/Tyr182) Antibody	Human Mouse Rat	WB
21584	PKC α Antibody	Human Mouse Rat	WB
11456	PKM2(phospho-Ser37) Antibody	Human	WB
21405	PP2A-alpha Antibody	Human Mouse Rat	WB
21409	Prospero homeobox protein 1 Antibody	Human Mouse	WB
21628	Rab5 Antibody	Human Mouse Rat	WB
21416	RE57611p Antibody	Fruit fly	WB
21682	S6 Ribosomal Protein (Ab-235/236) Antibody	Human Mouse Rat	WB
11580	S6 Ribosomal Protein (Phospho-Ser235/236) Antibody	Human Mouse Rat	WB
11579	Smad1 (Phospho-Ser206) Antibody	Human Mouse Rat	WB



21514 eIF4G (Ab-1232) Antibody



11577 Histone H3 (Phospho-Thr11) Antibody



11227 Met (Phospho-Tyr1234) Antibody



11221 mTOR (Phospho-Ser2448) Antibody