

Vertex[®]

The Peak of Pipetting

Filtered tips shown in actual size.



Vertex Macro Volume Tip Compatibility Chart

- Improved accuracy over serological pipettes when measuring large volumes of solution.
- Tip length facilitates access into bottles, cell culture flasks and tubes.
- Free of detectable RNase, DNase, DNA, and PCR inhibitors. Sterile versions are also pyrogen-free.

Macro Tips Description

Item #	Volume	Description	Packaging
4410-00	5 mL	Gilson-fit, Bulk, Natural	250/unit, 10 units/case
4420-00	5 mL	Eppendorf-fit, Bulk, Natural	
4411-00	5 mL	Gilson-fit, Racked, Natural	50 tips/unit, 10 units/case
4421-00	5 mL	Eppendorf-fit, Racked, Natural	
4411-S0	5 mL	Gilson-fit, Racked, Sterile, Natural	
4421-S0	5 mL	Eppendorf-fit, Racked, Sterile, Natural	
4411-SF	5 mL	Gilson-fit, Racked, Sterile, Filtered, Natural	
4421-SF	5 mL	Eppendorf-fit, Racked, Sterile, Filtered, Natural	
4520-01	10 mL	Bulk, Natural	200/unit, 5 units/case
4530-00	10 mL	Gilson-fit, Bulk, Natural	
4521-01	10 mL	Racked, Natural	25 tips/rack, 10 racks/case
4531-00	10 mL	Gilson-fit, Racked, Natural	
4521-S1	10 mL	Racked, Sterile, Natural	
4531-S0	10 mL	Gilson-fit, Racked, Sterile, Natural	
4521-SF	10 mL	Racked, Sterile, Filtered, Natural	
4531-SF	10 mL	Gilson-fit, Racked, Sterile, Filtered, Natural	



Macro Tips Compatibility

Item #	Eppendorf	Socorex	Gilson	Rainin	Accumax	Biohit	Finnpipette
4410-00			▲	▲			
4420-00	▲	▲					
4411-00			▲	▲			
4421-00	▲	▲					
4411-S0			▲	▲			
4421-S0	▲	▲					
4411-SF			▲	▲			
4421-SF	▲	▲					
4520-01		▲			▲	▲	▲
4530-00	▲		▲	▲			
4521-01		▲			▲	▲	▲
4531-00	▲		▲	▲			
4521-S1		▲			▲	▲	▲
4531-S0	▲		▲	▲			
4521-SF		▲			▲	▲	▲
4531-SF	▲		▲	▲			