

Background Information:

Nafuredin was isolated from Aspergillus niger which is marine strain. It inhibits NADH-fumarate reductase of anaerobic electron transport system in helminth mitochondria. Moreover, nafuredin exerts anthelmintic activity against Haemonchus contortus in in vivo trials with sheep.

Handling and Storage:

Store at -20°C. Do not dry up.

References:

1. J Antibiot (Tokyo). 2001 Mar;54(3):234-8.

2. Proc Natl Acad Sci U S A. 2001 Jan 2;98(1):60-2.