

Origin: Aspergillus sp.

CAS Registry Number: 483305-08-4

Appearance: White solid

Molecular Formula/ Weight:

C24H35NO4=401.55

Purity:>95.0% by HPLC

Solubility: Sol. in 1mg/ml ethanol, methanol,

DMSO

Background Information:

Aspergillin PZ, isoindole-alkaloid, was isolated from Aspergillus awamori. It induce morphological deformation of the conidia of Pyricularia oryzae.

Handling and Storage:

Store at -20°C.

References:

- 1. J Antibiot (Tokyo). 2002 Aug;55(8):693-5.
- 2. Tetrahedron. 2011 Dec 23;67(51):9837-9843.

Manufactured with Cortesy strain from The Kitasato Institute.