



Recombinat human Paraoxonase 1 (PON1)

ORDERING INFORMATION

Catalog Number: BML502

Lot Number:

Size: 10 µg

Expired Date: December 31, 2007

Formulation: PBS

Storage: -80°C

Specificity: human PON1 detected by specific monoclonal antibodies to PON1. Not analyzed enzyme activity

Application: Western blots, ELISA

Description

Recombinant human paraoxonase 1 (PON1) expressed in baculoviral insect expression system. The protein is suspended in phosphate buffered saline.

Purity

> 95% by SDS-PAGE.

Storage

-80°C. After initial defrost, aliquot products into individual tubes and refreeze at -80°C. Avoid repeated freeze/defrost cycles.

Stability

Protein is stable at least six month at -80°C without defrost.

Applications

Intended for use as a standard for western blots and ELISA.

References

(1) Kujiraoka et al., A sandwich enzyme-linked immunosorbent assay for human paraoxonase concentration. J Lipid Res, 2000;41:1355-1363.

FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

BML, Inc.

1361-1 Matoba, Kawagoe, Saitama 350-1101, Japan
TEL +81-49-232-0440/FAX +81-49-232-5480