

Anti-Medaka vitellogenin IgG, Horseradish Peroxidase conjugate

Code : MVD33-HRP

Lot :

Lot size : 0.5 mL

Background: Vitellogenin, the egg yolk precursor protein, is induced by estrogens.

On the other hand, environmental estrogens are known to disrupt endocrine system in animals. Recently, vitellogenin has proven to be an ideal marker for environmental estrogens.

Presentation : Supplied as a liquid in 0.01M Tris , 0.15M NaCl, 1% BSA, pH 7.5

Immunogen Source : Purified medaka vitellogenin

Applications : ELISA, Western blot

Normal dilution range : ELISA : 1:100 ~ 1:1000
Western blot : 1:100 ~ 1:5000

Storage : at -20°C for at least 12 months
Avoid freeze thaw cycles

Note : Since assay conditions vary, the optimum dilution should be determined for particular application.

Warning : For research use only

EnBioTec Laboratories, Co., Ltd.

Chiyoda-parion Bldg., 6th Floor, 2-3-16 Kandasuda-chou, Chiyoda-ku, Tokyo 101-0041, Japan

Tel:03-5297-7155 Fax:03-5297-0581 E-mail: support@enbiotec.co.jp