


IN3-VEC33 pFAST-G03



Name	Position	Name	Position
Left border (LB)	1-26	attR2	7163-7287
NOS terminator (NOST)	117-370	Intron	7355-7998
HPT	568-1590	attR2	8013-8137
NOS promoter (PNOS)	1605-1875	ccdB	8178-8483
35S terminator (T35S)	1950-2069	CmR	8825-9484
OLE1-GFP fusion protein (OLE1-GFP)	2353-3882	attR1	9593-9717
OLE1 promoter (POLE1)	3883-5326	35S promoter (P35S)	9829-10657
35S terminator (T35S)	5365-5484	Right border (RB)	10863-10887
attR1	5583-5707		
CmR	5816-6475		
ccdB	6817-7122		