

MOUSE ANTI-CA19.9  
MONOCLONAL ANTIBODY

CLONE:	3C11
QUANTITY:	0.1mg
SPECIFICITY:	Human CA19.9
IMMUNOGEN:	Purified CA19.9 from tumor cells
Type/Subtype:	IgA
APPLICATIONS:	Sandwich ELISA( This one clone is useful for both labeled antibody and capture antibody.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Salting -out fractionation by ammonium sulfate and ion exchange chromatographic purification)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8℃, Maintain at -70℃ for up to 12 months. Avoid repeated freeze/thaw cycle.

*For research use only; not for use as a diagnostic.*

Produced by:  
Japan Clinical Laboratories, Inc

Supplied by :  
Funakoshi Co., Ltd.