

MOUSE ANTI  $\alpha$  1-MICROGLOBULIN  
MONOCLONAL ANTIBODY

CLONE:	9F7
QUANTITY:	0.1mg
SPECIFICITY:	Human $\alpha$ 1-microglobulin
IMMUNOGEN:	Human $\alpha$ 1-microglobulin
Type/Subtype:	IgG1 $\kappa$
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a capture antibody in a human $\alpha$ 1-microglobulin ELISA in combination with biotinlated, monoclonal antibody 12D4.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

*For research use only; not for use as a diagnostic.*

Produced by:  
Japan Clinical Laboratories, Inc

Supplied by :  
Funakoshi Co., Ltd.