

MOUSE ANTI ANNEXIN V  
MONOCLONAL ANTIBODY

CLONE:	No.49
QUANTITY:	0.1mg
SPECIFICITY:	Human Annexin V
IMMUNOGEN:	Human Annexin V
Type/Subtype:	IgG2 b $\kappa$
APPLICATIONS:	Sandwich ELISA(This antibody can be used as a detection antibody in a human Annexin V ELISA combination with the monoclonal capture antibody No.23.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8℃, Maintain at -70℃ for up to 12 months. Avoid repeated freeze/thaw cycle.

*For research use only; not for use as a diagnostic.*

Produced by:  
Japan Clinical Laboratories,Inc

Supplied by :  
Funakoshi Co., Ltd.