

MOUSE ANTI FREE PSA
MONOCLONAL ANTIBODY

| | |
|-------------------|--|
| CLONE: | 2E2 |
| QUANTITY: | 0.1mg |
| SPECIFICITY: | Human Prostate specific antigen (PSA) of free form |
| IMMUNOGEN: | Purified antigen from human seminal fluid |
| Type/Subtype: | IgG1 κ |
| APPLICATIONS: | Free PSA sandwich ELISA (This antibody can be used as a capture antibody in a human Free PSA ELISA in combination with biotinlated, monoclonal antibody No.79. This clone specifically binds to free form of PSA and not binds to complex form.) |
| FORMAT: | Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography) |
| PRESENTATION: | Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2 |
| STORAGE/HANDLING: | 2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle. |

For research use only; not for use as a diagnostic.

Produced by:
Japan Clinical Laboratories, Inc

Supplied by :
Funakoshi Co., Ltd.