



Anti Human Organic Anion Transporter 2 (OAT2) Polyclonal Antibody


To eliminate the drug, xenobiotics, a variety of endogenous substances, and their metabolites out of the body, specific membrane proteins named transporters are required. There are two major pathways for the elimination, hepatic one through bile and renal one to urine. The transporter fall into various transport systems by the transportative substrate. In particular, organic ion transporter family is comprised of organic anion transport family (OAT), organic cation transport family (OCT), OCTN/carnitine transport family, and OAT are multispecific organic anion transporters, the substrates of which include a lot of both endogenous and exogenous anions .

Organic anion transporter 2 (OAT2) is predicted 12 putative membrane-spanning protein, and which is expressed predominantly in the liver and only weakly in the kidney. OAT2 mediates the Na^+ -independent transport of organic anions, such as salicylate, and also mediates the transport of prostaglandin E_2 (PGE_2), methotrexate, acetylsalicylate and PAH (ρ -aminohippurate).

This has been proved to be useful for immunohistochemistry .

| | |
|---|---|
| Package Size | 25 μg (250 μL / vial) |
| Format | Rabbit polyclonal antibody 0.1 mg/mL |
| Buffer | Block Ace as a stabilizer, containing 0.1%Proclin as bacteriostat |
| Storage | Store below -20°C Once thawed, store at 4°C . Repeated freeze-thaw cycles should be avoided. |
| Purification method | This antibody was purified from rabbit serum immunized with synthesized C-end peptide of human OAT2 by peptide affinity chromatography. |
| Working dilution for immunohistochemistry | 5 $\mu\text{g}/\text{mL}$ |
| HGNC Name | SLC722A7 (solute carrier family 22A7) |

*HGNC : Human Gene Nomenclature Committee



Human Kidney (frozen section)



Anti Human Organic Anion Transporter 2 (OAT2) Polyclonal Antibody

【Reference】

- 1 Sekine T., Cha S.H., Tsuda M., Apiwattanakul N., Nakajima N., Kanai Y., Endou H.:Identification of multispecific organic anion transporter 2 expressed predominantly in the liver. *FEBS Letters* 429 .179-182, 1998
- 2 Sekine T., Cha S.H., Kanai Y.and Endou H.:Molecular biology of multispecific organic anion transporter family (OAT family). *Clin.Exp.Nephrol*3.237-243, 1999
- 3 Sekine T., Cha S.H.and Endou H.:The multispecific organic anion transporter (OAT) family. *pflugers Arch-Eur.J.Physiol*440.337-350,2000

Manufacturer



Medicinal Chemistry Pharmaceutical Co., Ltd.

Kobe Research Institute

7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-945-7075 FAX:+81-78-306-0694

URL:<https://soyaku.co.jp/english/> tech-kobe@soyaku.co.jp

Previous manufacturer



Trans Genic Inc.