

## Anti human RXR gamma mouse monoclonal antibody

RXR gamma: Retinoid X Receptor gamma

Code No	PP-H3210-00	Application / Recommended Concentration In order to obtain the best results, optimal working dilutions should be determined by each individual user.	
Clone No.	H3210	Western Blot	1 ug/mL
Lot.	A-1	Non reducing Western Dist	
Concentration	1 mg/mL	Non reducing Western Blot	3 ug/mL
Volume	100 uL	ELISA	0.3 ug/mL
Ig Class	G2a	Immunoprecipitation	Decide by use
Description	Retinoid X receptor gamma (RXRg; NR2B3) is a member of orphan nuclear receptor. 9-cis retinoic acid can bind to RXR. RXRg is expressed brain, pituitary. Along with other members of the RXR family, RXRg plays roles in a variety of processes including embryonic patterning and organogenesis, cell proliferation and differentiation. RXRs commonly function as heterodimers with other members of the nuclear receptor superfamily. Two N-terminally different isoforms of RXRg with different expression patterns were found, RXRg1 and RXRg2.	Supershift Assay	Not yet tested
		Chromatin immunoprecipitation	Not yet tested
		Immunohistochemistry	Not yet tested
Nomenclature	NK2B3		
Genbank	U38480		

U38480		
Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human RXR gamma (2-137aa).	Storage	Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.
This antibody specifically recognizes human RXR gamma but does not recognize human RXR alpha and beta. Not yet tested in other species.	Reference	
Ammonium sulfate fractionation	-	
Physiological saline with 0.1% NaN3 as a preservative.	Notes	Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.
	Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human RXR gamma (2-137aa) . This antibody specifically recognizes human RXR gamma but does not recognize human RXR alpha and beta. Not yet tested in other species. Ammonium sulfate fractionation Physiological saline with 0.1% NaN3 as a	Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human RXR gamma (2-137aa) .       Storage         This antibody specifically recognizes human RXR gamma but does not recognize human RXR alpha and beta. Not yet tested in other species.       Reference         Ammonium sulfate fractionation       Notes

## FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited.

MADE IN JAPAN

Aug 28, 2006

## Distributor: Funakoshi co., Itd.

Funakoshi co.,Itd.
9-7 Hongo 2-chome, Bunkyo-ku, Tokyo 113-0033, Japan TEL: +81-3-5684-1620 FAX: +81-3-5684-1775 http://www.funakoshi.co.jp



Manufactured by **Perseus Proteomics Inc.** 4-7-6, Komaba, Meguro-ku, Tokyo 153-0041, Japan http://www.ppmx.com