

Anti human ERR beta mouse monoclonal antibody

ERR beta: Estrogen-Related Receptor beta

Code No	PP-H6707-00
Clone No.	H6707
Lot.	A-1
Concentration	1 mg/mL
Volume	100 uL
Ig Class	G2a
Description	Estrogen-related receptor beta (ERRb, ERR2; NR3B2) is a member of orphan nuclear receptor. ERRb is expressed in liver, brain, uterus, vagina and cervix. ERRb has a role in early placental development.

Application / Recommended Concentration

In order to obtain the best results, optimal working dilutions should be determined by each individual user.

Western Blot	Acceptant
western blot	1 ug/mL
Non reducing Western Blot	1 ug/mL
ELISA	0.2 ug/mL
Immunoprecipitation	Decide by use
Supershift Assay	Not yet tested
Chromatin immunoprecipitation	Not yet tested
Immunohistochemistry	Not yet tested

NR3B2	
X51417	=
Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human ERR beta(3-100 a.a.).	Stora
This antibody specifically recognizes human ERR beta and does not recognize human ERR alpha and gamma. Not yet tested in other species.	Refe
Ammonium sulfate fractionation	
Physiological saline with 0.1% NaN3 as a preservative.	- Note
	Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human ERR beta(3-100 a.a.). This antibody specifically recognizes human ERR beta and does not recognize human ERR alpha and gamma. Not yet tested in other species. Ammonium sulfate fractionation Physiological saline with 0.1% NaN3 as a

ige

Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.

rence

Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited. MADE IN JAPAN

Aug 22, 2006





9-7 Hongo 2-chome, Bunkyo-ku, TEL: +81-3-5684-1620 FAX: +81-3-5684-1775 http://www.funakoshi.co.jp



Manufactured by Perseus Proteomics Inc. 4-7-6, Komaba, Meguro-ku,

Tokyo 153-0041, Japan http://www.ppmx.com