

Anti human ERR gamma mouse monoclonal antibody

ERR gamma: Estrogen-Related Receptor gamma

Code No	PP-H6812-00
Clone No.	H6812
Lot.	A-1
Concentration	1 mg/mL
Volume	100 uL
lg Class	G2a
Description	Estrogen-related receptor gamma (ERRg, ESRRG; NR3B3) is a member of orphan nuclear receptor. ERRg is expressed in brain, kidney, testis, lung, adrenal gland, pancreas, placenta or bone marrow. ERRg was shown to be located in the critical region for Usher's syndrome.

Application / Recommended Concentration


In order to obtain the best results, optimal working dilutions should be determined by each individual user.

Western Blot	1 ug/mL
Non reducing Western Blot	Not yet tested
ELISA	1.7 ug/mL (A450=0.5)
Immunoprecipitation	Decide by use
Supershift Assay	Not yet tested

Chromatin immunoprecipitation Not yet tested

Immunohistochemistry

10µg/mL



Rat inner nuclear layer of the retina

Rat kidney distal convoluted tubule

with hybridoma of mouse myeloma cells (NS-1) and Storage Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated immunized with Baculovirus-expressed recombinant freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended. Reference Jun Qin, et al. Developmental Dynamics, 2007, 236: 810-820

> Notes Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts

of water during disposal.

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

preservative.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited. MADE IN JAPAN

Apr 17, 2008

Distributor: Funakoshi co., Itd.

Funakoshi co., Itd. 9-7 Hongo 2-chome, Bu Tokyo 113-0033, Japan 9-7 Hongo 2-chome, Bunkyo-ku, TEL: +81-3-5684-1620 FAX: +81-3-5684-1775 http://www.funakoshi.co.jp



Manufactured by Perseus Proteomics Inc. 4-7-6, Komaba, Meguro-ku, Tokyo 153-0041, Japan http://www.ppmx.com

human ERR gamma (2-100 aa) . Specificity This antibody specifically recognizes human ERR gamma and cross reacts with mouse and rat ERR gamma. This antibody does not recognize human ERR alpha and beta. Purification Ammonium sulfate fractionation Formulation Physiological saline with 0.1% NaN3 as a

Produced in BALB/c mouse ascites after inoculation

spleen cells derived from a BALB/c mouse

Nomenclature NR3B3

Genbank

Origin

AF058291